

Medical Image Analysis



॥ त्वं ज्ञानमयो विज्ञानमयोऽसि ॥

Angshuman Paul

Assistant Professor

Department of Computer Science & Engineering

Annotation-efficient DL for MedIA

Learning a New Concept: Human vs Computer

- Human learning
 - Usually requires only a few instances
 - May efficiently learn visual appearance from semantic description
- Machine Learning (especially deep learning)
 - May require thousands (if not millions) of instances
 - Difficult to train using semantic information for visual understanding

What is Annotation-efficient Machine Learning

- Learning new concepts from only a few or no labeled training data
- Mimicking the human cognitive ability
- Overcoming the curse of large annotation
- Small training set: Difficult to train deep models


Why Annotation-efficient Machine Learning in Medical Images?

- Learn new diseases from a few labeled examples
 - Rare disease diagnosis
- Working solutions with small labeled dataset
 - Reducing the need of manual annotations
 - Expert-dependent
 - Time consuming

Annotation-efficient Machine Learning




Few-shot
Learning



Zero-shot
Learning

Annotation-efficient Machine Learning for Chest X-ray Diagnosis: Challenges




Zero-shot Chest X-ray Diagnosis


Zero-Shot Learning for Medical Images

- Common diseases: Large pool of annotated training images
- Rare diseases: few or no annotated training images
- Use annotated images of common diseases
- Learn new diseases from no examples: use **auxiliary information**
 - Rare disease diagnosis


Zero-shot Learning




Zero-shot Learning: Training




Generalized Zero-shot Learning: Testing



Zero-shot Learning




Auxiliary Information for Radiology Diagnosis




Disease Signatures


- Signature Generator: Intelligent Word Embedding (IWE)¹




Zero-shot Diagnosis of Chest X-rays




Zero-shot Diagnosis of Chest X-rays




Zero-shot Diagnosis of Chest X-rays




Cross-Modality Semantic Embedding




Meaningful Semantic Embedding




Meaningful Semantic Embedding



Meaningful Semantic Embedding




Meaningful Semantic Embedding




$$L_{se} = \sum d - (p + q + r + s)$$


Meaningful Semantic Embedding



Meaningful Semantic Embedding




Meaningful Semantic Embedding



$$L_{sal} = \sum_i \left(\frac{1}{||v_i||^2 + \epsilon} \right)$$


Process Pipeline




Feature Extractor

- DenseNet-121²
- Trained using seen classes
- Feature vectors from the penultimate layer
- Feature vectors: likely to be **noisy** for unseen classes


Process Pipeline



Choice of Embedding Network



Embedding Network: Autoencoders³




- Input data x
- Hidden state z
- Reconstructed output \hat{x}


Embedding Network: Autoencoders

- Losses
 - Semantic embedding loss L_{se}
 - Semantic saliency loss L_{sal}
 - Reconstruction loss $L_{re} = ||x - \hat{x}||$
- Semantic Saliency Preserving Autoencoders (SSP-AE)

Dealing with Noisy Feature Vectors

- Ensemble⁴ of Autoencoders
- Several autoencoders trained in parallel
- Each autoencoder
 - Explores different feature subspaces
 - Semi-deterministic selection of subspaces
 - Trained with different bootstrap samples⁴











Seen & Unseen Classes


Class	
Pneumonia	← Seen Classes
Nodule	
Pneumothorax	
Consolidation	
Effusion	
Infiltration	
Edema	← Unseen Classes
Emphysema	
Cardiomegaly	

Visual Results

Unseen Class			Seen Class			
CXR Examples						
GT	Cardiomegaly	Infiltration Pneumothorax Emphysema	Cardiomegaly Edema	Nodule	Effusion Infiltration	Pneumothorax
D	Cardiomegaly	Emphysema	Nodule	Nodule	Effusion	Pneumonia

Can We Use More Auxiliary Information?


Trait-guided Multi-view Semantic Embedding for Zero-shot Chest X-ray Diagnosis



Visual Traits

- Location: Contains anatomical location information
 - lung, heart, etc.
- Position: Contains information of the pixel position corresponding to the abnormality
 - Upper portion of image, lower portion of image, etc.
- Opacity: high, low, medium
- Distribution: unilateral, bilateral
- Border sharpness: clear margin, indistinct margin
- Size (relative to lung volume)
- Aspect ratio

Semantic Spaces



Seen & Unseen Classes

Unseen Classes

Seen Classes

Class
Pneumonia
Nodule
Infiltration
Consolidation
Effusion
Pneumothorax
Edema
Emphysema
Cardiomegaly

Combination 1








Class
Pneumonia
Nodule
Infiltration
Consolidation
Effusion
Pneumothorax
Edema
Emphysema
Cardiomegaly

Combination 2

Class
Pneumonia
Nodule
Infiltration
Consolidation
Effusion
Pneumothorax
Edema
Emphysema
Cardiomegaly

Combination 3

Visual Results

Dataset	NIH-900		Open-i		PMC		
Image Examples							
Ground Truth	Cardiomegaly	Infiltration	Edema Cardiomegaly	Pneumonia	Consolidation	Nodule	Pneumothorax
Detected	Cardiomegaly (S)	Infiltration (U)	Edema (S)	Pneumonia (U)	Consolidation (S)	Nodule (U)	Nodule (U)

VAE for Zero-shot CXR Diagnosis








Figure 1: **Overview of our CXR-ML-GZSL model for learning visual representations of chest X-rays.** An overview of our network for chest X-ray images. It includes a trainable visual encoder and v - and d -dimensional visual and semantics spaces, respectively. For an input image x and its labels y , the network learns a visual representation guided by semantics extracted by BioBert. We perform end-to-end training of the visual encoder, visual mapping module, and the semantics mapping module, as indicated by the black dashed line.

Few-shot Chest X-ray Diagnosis Using Images from the Published Scientific Literature

Sample X-ray Images in Published Literature



Challenges in Few-shot Learning from Published Literature




- Learning from few images
- Artifacts, low-resolution


Solution: Use of Labeled and Unlabeled Images

- Labeled images from published literature: PubMed Central (PMC)
 - Initial training of few-shot learning model
- Unlabeled images: from NIH CXR dataset
 - High-resolution, less artifacts
- Re-training with pseudo labels for the NIH CXR dataset
 - Dealing with the problem of noisy labels


Process Pipeline




Process Pipeline: Initial Training of Classifier




The Loss Function




Loss Computed From Clusters



Loss Computed From Clusters



Loss Computed From Clusters




- Minimize every d_{ij} inside each cluster

$$L_{con} = \sum_{\forall c \in C} \left(\sum_{\forall x_i, x_j \in S_c} d_{ij} \right)$$

- Maximize r_{ab} for every pair of clusters


$$L_{sep} = - \sum_{a, b \in C} r_{ab}$$

Loss Computed From Clusters




- Minimize every d_{ij} inside each cluster (L_{con})
- Maximize r_{ab} for every pair of clusters (L_{sep})
- $L = L_{re} + \lambda_1 L_{con} + \lambda_2 L_{sep}$


Loss Computed From Clusters



Inference for Unlabeled Images: Pseudo Labels



Re-training



Inference













Image Results

Image Examples			
Ground Truth	Cardiomegaly	Edema	Edema
Detected	Cardiomegaly	Edema	Cardiomegaly


FSL Using an Ensemble of Subspaces




FSL Using an Ensemble of Subspaces

Images					
GT	Hernia	Fibrosis	Emphysema	Edema	Pneumonia
P	Hernia	Fibrosis	Emphysema	Cardiomegaly	Pneumonia

Few-shot Learning for Breast Cancer Detection




a) Binary classification on the ICIAR dataset by fine-tuning with ResNet50 V2 pre-trained architecture



b) Multi output model for primary and fine-grain classification trained on very few histopathology images from the BreakHis dataset. Model 1 has been used as the base model.

Lightweight CNN Model for Chest X-ray Diagnosis



~ 0.14M parameters and ~ 550 KB size

SI						
GT	Pleural Effusion	Cardiomegaly	Atelectasis	No Finding	Pneumothorax	Consolidation
D	Pleural Effusion	Cardiomegaly	Atelectasis	No Finding	Pneumothorax	Atelectasis

Single Image Super-resolution for Chest X-rays


Low Res Image




Our Model



Ground Truth



Federated Learning for Radiology Diagnosis



Dataset	NIH			CheXpert			VinBig		
Image Examples									
GT	Abnormal	Abnormal	Abnormal	Normal	Normal	Normal	Normal	Normal	Normal
Detected	Abnormal	Normal	Abnormal	Normal	Abnormal	Normal	Abnormal	Normal	Normal

Type	Method	F1 score			Accuracy		
		NIH	CX	VB	NIH	CX	VB
Local	NIH-L	0.93	0.89	0.64	0.93	0.81	0.68
	CX-L	0.90	0.90	0.70	0.90	0.84	0.75
	VB-L	0.80	0.87	0.88	0.83	0.79	0.90
Global	DN121-SA	0.90	0.89	0.70	0.91	0.82	0.79
	DN121-TL	0.83	0.88	0.65	0.85	0.81	0.78
	DN121-VL	0.88	0.88	0.69	0.89	0.81	0.80
	[6]	0.89	0.88	0.49	0.88	0.80	0.61
	Proposed	0.90	0.90	0.71	0.90	0.84	0.76